


ITS.

Implants
trauma



PROMET

Proximal Medial Tibia Locking Plate

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

WARNING: If there is no sufficient bone healing, wrong or incomplete postoperative care, plate might break.

All ITS plates are preformed anatomically as a matter of principle. If adjustment of the plate to the shape of the bone is required, this is possible by carefully bending gently in one direction once. Particular care is required when bending in the region of a plate hole, as deformation of the plate may lead to a failure of the locking mechanism. The plate must not be buckled or bent several times. This is particularly important in the case of titanium implants, to prevent material fatigue and subsequent failure. The method of bending is the conscious responsibility of the operating doctor; I.T.S. GmbH can accept no liability whatsoever for this.

Contents

I. Introduction

- P. 5 Preface
- P. 6 Screws
- P. 7 Properties
- P. 8 Pre-operative planning
- P. 8 Indications & Contraindications

2. Surgical technique

- P. 10 Time of operation
- P. 10 Pre-operative patient preparation
- P. 10 Access
- P. 10 Reduction
- P. 11 Placement of the screws
- P. 14 Postoperative treatment
- P. 14 Explantation

3. Information

- P. 16 Notes
- P. 17 Locking
- P. 17 Dotize®
- P. 18 Order list

Introduction



○ Preface

The Locking Proximal Medial Tibia Plate is an osteosynthesis system for various proximal tibia fractures, and can be used for complex C1 to C3 fractures in combination with the Proximal Lateral Tibia Plate.

The special feature of this implant is the free choice of screw placement. The user is able to set any desired screw in any hole, either locking or non-locking screw (except oblong hole).

The free choice of screw angulation ($\pm 15^\circ$, see page 17) as well as the possibility to set several subchondral screws provides advantages especially in the case of complex fractures.



○ Screws

3735I-XX-N Cortical Screw, Locking, D=3.5mm, SH

61273-220 Spiral Drill, D=2.7mm, L=220mm, AO Connector


56252 Screwdriver, WS 2.5,
self-holding sleeve



3235I-XX Cortical Screw, D=3.5mm

61273-220 Spiral Drill, D=2.7mm, L=220mm, AO Connector

56252 Screwdriver, WS 2.5,
self-holding sleeve



37422-XX-N Cancellous Screw, Locking, D=4.2mm, SH

61253-220 Spiral Drill, D=2.5mm, L=220mm, AO Connector

56252 Screwdriver, WS 2.5,
self-holding sleeve



35164-I50 Guide Wire, Steel, D=1.6mm,
L=150mm, TR, w. thread




○ Properties

Properties of the material:


- Plate material: Titanium
- Material of screws: TiAl6V4 ELI
- Easier removal of the implant after the fracture has healed
- Improved fatigue strength of the implant
- Reduced risk of cold welding
- Reduced risk of inflammation and allergy

Properties of the implant:

- Multi-directional Locking
- Anatomical plate design
- 5 proximal plate holes for peri-articular fixation
- Oblong hole for optimal positioning and alignment of the tibia length
- Pointed distal plate end for percutaneous insertion



○ Pre-operative planning



○ Indications, Contraindications

Indications:

- Buttress metaphyseal fractures of the medial tibia plateau
- Split fractures of the medial tibia plateau (with possible depressions)
- Stabilization of the proximal quarter of the tibia (medial and lateral)
- Complex CI-C3 fractures can be treated in combination with the PTL - Proximal Lateral Tibia Plate

Contraindications:

- Existing infections in the fracture zone and operation area
- Common situations that do not allow osteosynthesis
- With advanced osteoporosis
- In cases of skin and soft tissue problems that prevent a tension-free skin closure
- Obesity
- Lack of patient compliance

2.

◦ Time of operation

- ♦ Immediately after trauma or delayed

◦ Pre-operative patient preparation

- ♦ General anaesthesia, local anaesthesia or combination can be used
- ♦ The patient is in the supine position with the leg slightly in flexion on a pedestal
- ♦ Application of a tourniquet


◦ Access

Posteromedial incision:

- ♦ Straight or slightly curved incision between the hamstrings (pes anserinus) and gastrocnemius muscle, running from the medial epicondyle towards the posteromedial edge of the tibia.

◦ Reduction

- ♦ Temporary fixation of the plate using guide wires or optionally available Temporary Plate Holder (58164-I50).
- ♦ Anatomical reduction of the articular block and fracture segments to the plate (varus/valgus, ante-/retroversion)
- ♦ Subsequent control under fluoroscopy



○ Placement of the screws

For optimal alignment of the length of reduction, we recommend to first fill the oblong hole.


With the spiral drill, D=2.7mm, L=220mm, AO Connector (**61273-220**), drill through the drill guide, D=2.7/2.0mm (**62202**) into the oblong hole.

Determine appropriate length using the depth gauge, solid small fragment screws (**59022**). Insert the D=3.5mm cortical screw (**3235I-XX**) with the screwdriver, WS 2.5, self-holding sleeve (**56252**).



Then, using the spiral drill, D=2.5mm, L=220mm, AO Connector (**61253-220**), drill through the drill guide, D=2.7/2.0mm (**62202**) into a proximal plate hole.

Determine appropriate length using the depth gauge, solid small fragment screws (**59022**). Insert the D=4.2mm locking cancellous screw (**37422-XX-N**) with the screwdriver, WS 2.5, self-holding sleeve (**56252**).



Using the spiral drill, D=2.7mm, L=220mm, AO Connector (**61273-220**) drill through the drill guide, D=2.7/2.0mm (**62202**) into a shaft plate hole.

Determine appropriate length using the depth gauge, solid small fragment screws (**59022**). Insert the D=3.5mm cortical screw, optionally locking (**32351-XX / 37351-XX-N**) with the screwdriver, WS 2.5, self-holding sleeve (**56252**).



The remaining plate holes are then filled, with either locking or non-locking screws. Subsequent control of plate position under fluoroscopy.



◦ Postoperative treatment

- Keep leg raised for 2 to 5 days and take decongestant actions
- After reduction of swelling, beginning of the passive mobilization (CPM splint)
- No weight bearing in the treatment of articular fractures for a minimum of 10 to 12 weeks
- When a locking screw connection has been used, it is necessary to be aware that the diagnosis of a non-union may be very delayed.

◦ Explantation

If desired by the patient, the implant can be removed. Removal should be performed at the earliest 1 1/2 years later or after radiographic verification of the healed bone.

The problem of cold welding was resolved by using a special surface treatment (for further information see page 17).

3.


Notes

[illegible]

○ Locking

Locking works because:

- Screw material (TiAlV) is slightly harder than plate material (Titanium Grade 2)
- Screw head **forms** thread into the plate (no cutting)



Benefits:

- $\pm 15^\circ$ and Locking
- No pre threading
- No cold welding
- No debris
- You can re-set the screw up to 3 times

○ Dotize®


Chemical process - anodization in a strong alkaline solution*

Type III anodization

- Layer thickness 60-200nm
- + Different colors
- Implant surface remains sensitive to:
Chipping
Peeling
Discoloration

Dotize Type II anodization

- Layer thickness 2000-10 000nm
- + Film becomes an interstitial part of the titanium
- No visible cosmetic effect



Anodization Type II leads to following benefits*

- Oxygen and silicon absorbing conversion layer
- Decrease in protein adsorption
- Closing of micro pores and micro cracks
- Reduced risk of inflammation and allergy
- Hardened titanium surface
- Reduced tendency of cold welding of titanium implants
- Increased fatigue resistance of implants
- Improved wear and friction characteristics

Order list

Proximal Medial Tibia Plate, 4-Hole
Proximal Medial Tibia Plate, 7-Hole

21325-4
21325-7



Cortical Screw, D=3.5mm, L=24mm
Cortical Screw, D=3.5mm, L=28mm
Cortical Screw, D=3.5mm, L=32mm
Cortical Screw, D=3.5mm, L=36mm
Cortical Screw, D=3.5mm, L=40mm
Cortical Screw, D=3.5mm, L=42mm
Cortical Screw, D=3.5mm, L=44mm
Cortical Screw, D=3.5mm, L=46mm
Cortical Screw, D=3.5mm, L=48mm
Cortical Screw, D=3.5mm, L=50mm
Cortical Screw, D=3.5mm, L=55mm
Cortical Screw, D=3.5mm, L=60mm
Cortical Screw, D=3.5mm, L=65mm
Cortical Screw, D=3.5mm, L=70mm
Cortical Screw, D=3.5mm, L=75mm
Cortical Screw, D=3.5mm, L=80mm
Cortical Screw, D=3.5mm, L=85mm
Cortical Screw, D=3.5mm, L=90mm

32351-24
32351-28
32351-32
32351-36
32351-40
32351-42
32351-44
32351-46
32351-48
32351-50
32351-55
32351-60
32351-65
32351-70
32351-75
32351-80
32351-85
32351-90



Cortical Screw, Locking, D=3.5mm, L=24mm, SH
Cortical Screw, Locking, D=3.5mm, L=26mm, SH
Cortical Screw, Locking, D=3.5mm, L=28mm, SH
Cortical Screw, Locking, D=3.5mm, L=30mm, SH
Cortical Screw, Locking, D=3.5mm, L=32mm, SH
Cortical Screw, Locking, D=3.5mm, L=34mm, SH
Cortical Screw, Locking, D=3.5mm, L=36mm, SH
Cortical Screw, Locking, D=3.5mm, L=38mm, SH
Cortical Screw, Locking, D=3.5mm, L=40mm, SH
Cortical Screw, Locking, D=3.5mm, L=42mm, SH
Cortical Screw, Locking, D=3.5mm, L=44mm, SH
Cortical Screw, Locking, D=3.5mm, L=46mm, SH
Cortical Screw, Locking, D=3.5mm, L=48mm, SH
Cortical Screw, Locking, D=3.5mm, L=50mm, SH
Cortical Screw, Locking, D=3.5mm, L=55mm, SH
Cortical Screw, Locking, D=3.5mm, L=60mm, SH

37351-24-N
37351-26-N
37351-28-N
37351-30-N
37351-32-N
37351-34-N
37351-36-N
37351-38-N
37351-40-N
37351-42-N
37351-44-N
37351-46-N
37351-48-N
37351-50-N
37351-55-N
37351-60-N



Cancellous Screw, Locking, D=4.2mm, L=18mm, SH
Cancellous Screw, Locking, D=4.2mm, L=20mm, SH
Cancellous Screw, Locking, D=4.2mm, L=22mm, SH
Cancellous Screw, Locking, D=4.2mm, L=24mm, SH
Cancellous Screw, Locking, D=4.2mm, L=26mm, SH
Cancellous Screw, Locking, D=4.2mm, L=28mm, SH
Cancellous Screw, Locking, D=4.2mm, L=30mm, SH
Cancellous Screw, Locking, D=4.2mm, L=32mm, SH
Cancellous Screw, Locking, D=4.2mm, L=34mm, SH
Cancellous Screw, Locking, D=4.2mm, L=36mm, SH

37422-18-N
37422-20-N
37422-22-N
37422-24-N
37422-26-N
37422-28-N
37422-30-N
37422-32-N
37422-34-N
37422-36-N



Cancellous Screw, Locking, D=4.2mm, L=38mm, SH	37422-38-N
Cancellous Screw, Locking, D=4.2mm, L=40mm, SH	37422-40-N
Cancellous Screw, Locking, D=4.2mm, L=42mm, SH	37422-42-N
Cancellous Screw, Locking, D=4.2mm, L=44mm, SH	37422-44-N
Cancellous Screw, Locking, D=4.2mm, L=46mm, SH	37422-46-N
Cancellous Screw, Locking, D=4.2mm, L=48mm, SH	37422-48-N
Cancellous Screw, Locking, D=4.2mm, L=50mm, SH	37422-50-N
Cancellous Screw, Locking, D=4.2mm, L=55mm, SH	37422-55-N
Cancellous Screw, Locking, D=4.2mm, L=60mm, SH	37422-60-N
Cancellous Screw, Locking, D=4.2mm, L=65mm, SH	37422-65-N
Cancellous Screw, Locking, D=4.2mm, L=70mm, SH	37422-70-N
Cancellous Screw, Locking, D=4.2mm, L=75mm, SH	37422-75-N
Cancellous Screw, Locking, D=4.2mm, L=80mm, SH	37422-80-N
Cancellous Screw, Locking, D=4.2mm, L=85mm, SH	37422-85-N
Cancellous Screw, Locking, D=4.2mm, L=90mm, SH	37422-90-N

Screwdriver, WS 2.5, self-holding sleeve

56252



Depth Gauge, Solid Small Fragment Screws

59022



Drill Guide, D=2.0/2.7mm

62202



Spiral Drill, D=2.5mm, L=220mm, AO Connector

61253-220

Spiral Drill, D=2.7mm, L=220mm, AO Connector

61273-220



Guide Wire, Steel, D=1.6mm, L=150mm, TR, w. Thread

35164-150



Sterilization Tray, Proximal Medial Tibia Plate

50280

Optionally (on request)

Temporary Plate Holder

58164-150



For detailed cleaning and sterilization instructions, please refer to package insert.

* All implants are available sterile-packed optionally.
Add „-S“ to article number for sterile-packed implants (ex. 37304-12-S; 21031-3-S)
Delivery times, prices & minimum quantity vary from standard.



I.T.S. USA
1778 Park Avenue N, Suite 200
Maitland, FL 32751

Tel.: 877 - 971 - 8054
Fax: 877 - 971 - 8056
info@its-implantusa.com
www.its-implantusa.com



Order No. PROMET-OP-0717-USA
Edition: July/2017

© ITS. GmbH Graz/Austria 2017.
Subject to technical alterations, errors and misprints excepted.