


ITS.

Implants
trauma



CAL

Calcaneus Locking Plate

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

WARNING: If there is no sufficient bone healing, wrong or incomplete postoperative care, plate might break.

All ITS plates are preformed anatomically as a matter of principle. If adjustment of the plate to the shape of the bone is required, this is possible by carefully bending gently in one direction once. Particular care is required when bending in the region of a plate hole, as deformation of the plate may lead to a failure of the locking mechanism. The plate must not be buckled or bent several times. This is particularly important in the case of titanium implants, to prevent material fatigue and subsequent failure. The method of bending is the conscious responsibility of the operating doctor; I.T.S. GmbH can accept no liability whatsoever for this.

Contents

I. Introduction

- P. 5 Preface
- P. 6 Screws
- P. 7 Properties
- P. 8 Indications & Contraindications
- P. 8 Time of operation

2. Surgical Technique

- P. 10 Pre-operative patient preparation
- P. 10 Calcaneus Distractor & Heel Press
- P. 10 Diagnosis
- P. 11 Access
- P. 11 Reduction
- P. 12 Fixation
- P. 13 Postoperative treatment
- P. 13 Explantation

3. Information

- P. 15 Locking
- P. 15 Dotize®
- P. 16 Order list

Introduction



○ Preface

The Locking Calcaneus Plate is a proven osteosynthesis system for various calcaneus fractures.

The special feature of this implant is the free choice of screw placement. The user is able to set any desired screw in any hole (either using a locking or non-locking screw).

The free choice of screw angulation ($\pm 15^\circ$, see page 15) provides an advantage in fracture treatment, especially in the case of complex fractures.





○ Screws

3235I-XX Cortical Screw, D=3.5mm

61273-I00 Spiral Drill, D=2.7mm, L=100mm, AO Connector



56252 Screwdriver, WS 2.5,
with self-holding sleeve



37422-XX-N Cancellous Screw, locking, D=4.2mm, SH

61253-I10 Spiral Drill, D=2.5mm, L=110mm, AO Connector

56252 Screwdriver, WS 2.5,
with self-holding sleeve




3735I-XX Cortical Screw, Locking, D=3.5mm, SH

61273-I00 Spiral Drill, D=2.7mm, L=100mm, AO Connector

56252 Screwdriver, WS 2.5,
with self-holding sleeve



OPTIONAL
(ON REQUEST)




○ Properties

Properties of the material:

- ◆ Plate material: Titanium
- ◆ Material of screws: TiAl6V4 ELI
- ◆ Easier removal of the implant after the fracture has healed
- ◆ Improved fatigue strength of the implant
- ◆ Reduced risk of cold welding
- ◆ Reduced risk of inflammation and allergy

Properties of the implant:

- ◆ Multi-directional locking
- ◆ Anatomical plate design
- ◆ Short and long version



◦ Indications, Contraindications & Time of operation

Indications:

- ♦ Complex fractures of the calcaneus
- ♦ All intra-articular fractures with relevant joint distortion and a comminution zone in which a semi-operative procedure (screws, drill wires) does not raise expectations of exact repositioning

Contraindications:

- ♦ Existing infections in the fracture zone and operation area
- ♦ Common situations that do not allow osteosynthesis
- ♦ Obesity
- ♦ Lack of patient compliance

Time of operation:

- ♦ After regression of the swelling

2.

◦ Pre-operative patient preparation

- ◊ Supine position or lateral position
- ◊ Blood arrest


◦ Diagnosis

Standard X-ray of the calcaneus, axial and coronet CT with reconstructions.

○ Access

Expanded lateral approach:

- Subperiosteal single layered lifting of a lateral skin-soft tissue flap
- Hold away the flap by using bent guide wires






○ Reduction

- Open reduction under view by mean of Schanz screw, Steinmann nails, guide wires
- Padding of comminuted zones with bone replacement

◦ Fixation

- ♦ Temporarily by guide wires, eventually by cancellous screws
- ♦ Forming the locking plate for calcaneus
- ♦ Fixation with the cancellous screws in the main fragments
- ♦ Finishing with locking screws



○ Postoperative treatment

- Plastex cast of the lower leg for 2 weeks until healing
- Physical therapy
- Mobilization by crutches
- Relief of the strain for 8-12 weeks

○ Explantation

If desired by the patient, the implant can be removed.

Removal should be performed at the earliest 1 1/2 years later or after radiographic verification of the healed bone.

The problem of cold welding was resolved by using a special surface treatment (for further information see page 15).


Information

3.

○ Locking

Locking works because:

- Screw material (TiAlV) is slightly harder than plate material (Titanium Grade 2)
- Screw head **forms** thread into the plate (no cutting)



Benefits:

- $\pm 15^\circ$ and Locking
- No pre threading
- No cold welding
- No debris
- You can re-set the screw up to 3 times

○ Dotize®


Chemical process - anodization in a strong alkaline solution*

Type III anodization

- Layer thickness 60-200nm
- + Different colors
- Implant surface remains sensitive to:
Chipping
Peeling
Discoloration

Dotize Type II anodization

- Layer thickness 2000-10 000nm
- + Film becomes an interstitial part of the titanium
- No visible cosmetic effect



Anodization Type II leads to following benefits*

- Oxygen and silicon absorbing conversion layer
- Decrease in protein adsorption
- Closing of micro pores and micro cracks
- Reduced risk of inflammation and allergy
- Hardened titanium surface
- Reduced tendency of cold welding of titanium implants
- Increased fatigue resistance of implants
- Improved wear and friction characteristics

Order list

Calcaneus Plate, Short
Calcaneus Plate, Long

21061-51
21061-60



| | |
|--|------------|
| Cancellous Screw, Locking, D=4.2mm, L=26mm, SH | 37422-26-N |
| Cancellous Screw, Locking, D=4.2mm, L=28mm, SH | 37422-28-N |
| Cancellous Screw, Locking, D=4.2mm, L=30mm, SH | 37422-30-N |
| Cancellous Screw, Locking, D=4.2mm, L=32mm, SH | 37422-32-N |
| Cancellous Screw, Locking, D=4.2mm, L=34mm, SH | 37422-34-N |
| Cancellous Screw, Locking, D=4.2mm, L=36mm, SH | 37422-36-N |
| Cancellous Screw, Locking, D=4.2mm, L=38mm, SH | 37422-38-N |
| Cancellous Screw, Locking, D=4.2mm, L=40mm, SH | 37422-40-N |
| Cancellous Screw, Locking, D=4.2mm, L=42mm, SH | 37422-42-N |
| Cancellous Screw, Locking, D=4.2mm, L=44mm, SH | 37422-44-N |
| Cancellous Screw, Locking, D=4.2mm, L=46mm, SH | 37422-46-N |
| Cancellous Screw, Locking, D=4.2mm, L=48mm, SH | 37422-48-N |
| Cancellous Screw, Locking, D=4.2mm, L=50mm, SH | 37422-50-N |



| | |
|---------------------------------|----------|
| Cortical Screw, D=3.5mm, L=26mm | 32351-26 |
| Cortical Screw, D=3.5mm, L=28mm | 32351-28 |
| Cortical Screw, D=3.5mm, L=30mm | 32351-30 |
| Cortical Screw, D=3.5mm, L=32mm | 32351-32 |
| Cortical Screw, D=3.5mm, L=34mm | 32351-34 |
| Cortical Screw, D=3.5mm, L=36mm | 32351-36 |
| Cortical Screw, D=3.5mm, L=38mm | 32351-38 |
| Cortical Screw, D=3.5mm, L=40mm | 32351-40 |
| Cortical Screw, D=3.5mm, L=42mm | 32351-42 |
| Cortical Screw, D=3.5mm, L=44mm | 32351-44 |
| Cortical Screw, D=3.5mm, L=46mm | 32351-46 |
| Cortical Screw, D=3.5mm, L=48mm | 32351-48 |
| Cortical Screw, D=3.5mm, L=50mm | 32351-50 |



Screwdriver, WS 2.5, with self-holding sleeve 56252



Depth Gauge, Solid Small Fragment Screws 59022



Drill Guide, D=2.0/2.7mm 62202



Spiral Drill, D=2.5mm, L=110mm, AO Connector 61253-110
Spiral Drill, D=2.7mm, L=100mm, AO Connector 61273-100



Sterilization Tray, Calcaneus Plate 50171

For detailed cleaning and sterilization instructions, please refer to package insert.

Special sizes & instruments optional on request *

| | |
|--|----------|
| Cortical Screw, Locking, D=3.5mm, L=26mm, SH | 32351-26 |
| Cortical Screw, Locking, D=3.5mm, L=28mm, SH | 32351-28 |
| Cortical Screw, Locking, D=3.5mm, L=30mm, SH | 32351-30 |
| Cortical Screw, Locking, D=3.5mm, L=32mm, SH | 32351-32 |
| Cortical Screw, Locking, D=3.5mm, L=34mm, SH | 32351-34 |
| Cortical Screw, Locking, D=3.5mm, L=36mm, SH | 32351-36 |
| Cortical Screw, Locking, D=3.5mm, L=38mm, SH | 32351-38 |
| Cortical Screw, Locking, D=3.5mm, L=40mm, SH | 32351-40 |
| Cortical Screw, Locking, D=3.5mm, L=42mm, SH | 32351-42 |
| Cortical Screw, Locking, D=3.5mm, L=44mm, SH | 32351-44 |
| Cortical Screw, Locking, D=3.5mm, L=46mm, SH | 32351-46 |
| Cortical Screw, Locking, D=3.5mm, L=48mm, SH | 32351-48 |
| Cortical Screw, Locking, D=3.5mm, L=50mm, SH | 32351-50 |



Tray Setting



Notes

This image shows a full page of blank handwriting practice paper. It features approximately 20 evenly spaced horizontal red lines across the entire page, providing a guide for letter height and placement. The background is plain white, and there are no margins, text, or other markings present.



I.T.S. USA
1778 Park Avenue N, Suite 200
Maitland, FL 32751

Tel.: 877 - 971 - 8054
Fax: 877 - 971 - 8056
info@its-implantusa.com
www.its-implantusa.com



Order No. CAL-OP-0218-USA
Edition: February/2018

© ITS. GmbH Graz/Austria 2018.
Subject to technical alterations, errors and misprints excepted.