




	Spongiosa- schraube	Spongiosa- schraube	Kortikalis- schraube	Stabilisie- rungs- schraube	Spongiosa Stabilisie- rungs- schraube	Kortikalis Stabilisie- rungs- schraube	Kortikalis- schraube	Kortikalis- schraube	Spongiosa- schraube	Spongiosa- schraube	Kortikalis- schraube	Kortikalis- schraube	Kortikalis- schraube	Spongiosa- schraube	Kortikalis- schraube	Spongiosa- schraube
	winkel- stabil Ø 1.5mm	winkel- stabil Ø 2.3mm	Ø 1.8mm	winkel- stabil Ø 2.4mm	winkel- stabil Ø 3.0mm	winkel- stabil Ø 3.0mm	Ø 2.7mm	winkel- stabil Ø 3.5mm	winkel- stabil Ø 3.5mm	winkel- stabil Ø 4.2mm	Ø 3.5mm	Ø 3.5mm	winkel- stabil Ø 4.5mm	winkel- stabil Ø 5.9mm	Ø 4.5mm	Ø 5.9mm
	3715I-xx	3725I-xx	3215I-xx	3724I-xx	37303-xx	37304-xx	3227I-xx	3735I-xx-N	37352-xx-N	37422-xx-N	3235I-xx	32352-xx	37455-xx	37592-xx	32455-xx	3059I-xx
SPIRALBOHRER	Ø 1.1mm	Ø 1.8mm	Ø 1.5mm	Ø 1.8mm	Ø 2.0mm	Ø 2.4mm	Ø 2.0mm	Ø 2.7mm	Ø 2.0mm	Ø 2.5mm	Ø 2.7mm	Ø 2.7mm	Ø 3.2mm	Ø 3.5mm	Ø 3.2mm	Ø 3.5mm
Hand Locking Plates System 																
Ulna Osteotomy Locking Plate II 																
Distal Ulna Locking Plate 																
PROlock Radius Locking Plate II 																
PROlock Radius Locking Plate PHOENIX 																
Tubular Locking Plates System 																
Hallux Osteotomy Locking Plate 																
Foot Locking Plates System 																
Clavicle Locking Plate System NARROW 																
Clavicle Locking Plate System 																
Proximal Flanges Humeral Locking Plate PHOENIX 																
Proximal Standard Humeral Locking Plate PHOENIX 																
Distal Humeral Locking Plates 																
Olecranon Locking Plate System 																
Straight Compression Locking Plate 																
Straight Locking Plate System* * ausgenommen Straight Locking Plate 3.5/4.5mm 																
Universal Forearm Locking Plates System 																
Pelvic Small Fragment Reconstruction System PHOENIX 																
PROximal MEDial Tibia Locking Plate 																
Pilon Locking Plates System 																
Distal Anterolateral Tibia Locking Plate 																
Fibula Locking Plate 																
CALcaneus Locking Plate 																
Straight Locking Plate 3.5/4.5mm 																
Pelvic Large Fragment Reconstruction System PHOENIX 																
Distal Femur Locking Plate 																
Proximal Lateral Tibia Locking Plate 																
SCHRAUBENDREHER	Torque T5	Torque T9	Sechskant SW 2.5			Sechskant SW 3.5										