

		Spongiosa- schraube	Spongiosa- schraube	Kortikalis- schraube	Stabilisie- rungs- schraube	Spongiosa Stabilisie- rungs- schraube	Kortikalis Stabilisie- rungs- schraube	Kortikalis- schraube	Kortikalis- schraube	Spongiosa- schraube	Spongiosa- schraube	Kortikalis- schraube	Kortikalis- schraube	Kortikalis- schraube	Spongiosa- schraube	Kortikalis- schraube	Spongiosa- schraube
		winkel- stabil Ø 1.5mm	winkel- stabil Ø 2.3mm	Ø 1.8mm	winkel- stabil Ø 2.4mm	winkel- stabil Ø 3.0mm	winkel- stabil Ø 3.0mm	Ø 2.7mm	winkel- stabil Ø 3.5mm	winkel- stabil Ø 3.5mm	winkel- stabil Ø 4.2mm	Ø 3.5mm	Ø 3.5mm	winkel- stabil Ø 4.5mm	winkel- stabil Ø 5.9mm	Ø 4.5mm	Ø 5.9mm
		3715l-xx	3725l-xx	3215l-xx	3724l-xx	37303-xx	37304-xx	3227l-xx	3735l-xx-N	37352-xx-N	37422-xx-N	3235l-xx	32352-xx	37455-xx	37592-xx	32455-xx	3059l-xx
SPIRALBOHRER		Ø 1.1mm	Ø 1.8mm	Ø 1.5mm	Ø 1.8mm	Ø 2.0mm	Ø 2.4mm	Ø 2.0mm	Ø 2.7mm	Ø 2.0mm	Ø 2.5mm	Ø 2.7mm	Ø 2.7mm	Ø 3.2mm	Ø 3.5mm	Ø 3.2mm	Ø 3.5mm
Hand Locking Plates System																	
Ulna Osteotomy Locking Plate II																	
Distal Ulna Locking Plate																	
PROlock Radius Locking Plate II																	
PROlock Radius Locking Plate PHOENIX																	
Tubular Locking Plates System																	
Hallux Osteotomy Locking Plate																	
Foot Locking Plates System																	
Anterior Clavicle Locking Plates System																	
Clavicle Locking Plate System NARROW																	
Clavicle Locking Plate System																	
Proximal Flanges Humeral Locking Plate PHOENIX																	
Proximal Standard Humeral Locking Plate PHOENIX																	
Distal Humeral Locking Plates																	
Olecranon Locking Plate System																	
Straight Compression Locking Plate																	
Straight Locking Plate System*																	
* ausgenommen Straight Locking Plate 3.5/4.5mm																	
Universal Forearm Locking Plates System																	
Pelvic Small Fragment Reconstruction System PHOENIX																	
PROximal MEDial Tibia Locking Plate																	
Pilon Locking Plates System																	
Distal Anterolateral Tibia Locking Plate																	
Fibula Locking Plate																	
CALcaneus Locking Plate																	
CALcaneus mi Locking Plate																	
Straight Locking Plate 3.5/4.5mm																	
Pelvic Large Fragment Reconstruction System PHOENIX																	
Distal Femur Locking Plate																	
Proximal Lateral Tibia Locking Plate																	
SCHRAUBENDREHER		Torque T5															