


ITS.

Implants for Trauma Surgery



PROlock Radius Locking Plate

PHOENIX

STYLOID DESIGN

PROlock Radius Locking System PHOENIX



Screws:

- 3.0mm Locking Cancellous Screw (37303-xx)
- 2.7mm Cortical Screw (3227I-xx)


(Optional):

- 3.0mm Locking Cortical Screw (37304-xx)
- 2.4mm Locking Stabilization Screw (3724I-xx)



WATERSHED LINE DESIGN

PROlock Radius Locking System II



Screw Trajectories

PROlock Radius Locking System PHOENIX

26mm Plate								
	1	2	3	4	5	6	7	8
Proximal	0°	0°	0°	2°	0°	0°	0°	0°
Distal	15°	12°	0° (5°)	0°	12°	13°	4° (8°)	3°
Ulnar	0°	0°	0°	0°	0°	0°	0°	5°
Radial	8°	7°	0°	0°	7°	0°	0°	0°



Disclaimer:

The intended users are limited to medical personnel with appropriate product training by the medical product consultants or knowledge of the surgical procedure to be applied. The medical staff must ensure that the use of I.T.S. GmbH medical devices is appropriate, taking into account the medical condition and medical history of the patient. Prior to product use, medical personnel must refer to complete information on product label and in IFU, including, but not limited to, indications, contraindications, warnings and preventative measures, and cleaning and sterilization instructions. Product availability is dependent on country registrations and clearances. For more information, please visit www.its-implant.com or contact us at office@its-implant.com. All information herein is the intellectual property of I.T.S. GmbH.

Order Nr. PRLPHX-IP-0823-EN
Edition: August/2023



HEADQUARTER

I.T.S. GmbH

Autal 28, 8301 Lassnitzhöhe, Austria

Tel.: +43 (0) 316/ 211 21 0

office@its-implant.com

www.its-implant.com

I.T.S. Deutschland GmbH

Thurn-und-Taxis-Platz 6, 60313 Frankfurt

Tel.: +49 69 25 73 75 322

germany@its-implant.com

www.its-implant.com

I.T.S. USA

1778 Park Avenue N - Suite 200

Tel.: 407-971-8054

info@its-implantusa.com

www.its-implant.com

ITS.

Implants for Trauma Surgery



www.its-implant.com



ITS. GmbH



ITS. Implants



[its.implants](https://www.instagram.com/its.implants)