

Cancellous Screw	Cancellous Screw	Cortical Screw	Stabilization Screw	Cancellous Stabilization Screw	Cortical Stabilization Screw	Cortical Screw	Cortical Screw	Cancellous Screw	Cancellous Screw	Cortical Screw	Cortical Screw	Cortical Screw	Cancellous Screw	Cortical Screw	Cancellous Screw
locking Ø 1.5mm	locking Ø 2.3mm	Ø 1.8mm	locking Ø 2.4mm	locking Ø 3.0mm	locking Ø 3.0mm	Ø 2.7mm	locking Ø 3.5mm	locking Ø 3.5mm	locking Ø 4.2mm	Ø 3.5mm	Ø 3.5mm	locking Ø 4.5mm	locking Ø 5.9mm	Ø 4.5mm	Ø 5.9mm
3715I-xx	3725I-xx	3215I-xx	3724I-xx	37303-xx	37304-xx	3227I-xx	3735I-xx-N	37352-xx-N	37422-xx-N	3235I-xx	32352-xx	37455-xx	37592-xx	32455-xx	3059I-xx
Ø 1.1mm	Ø 1.8mm	Ø 1.5mm	Ø 1.8mm	Ø 2.0mm	Ø 2.4mm	Ø 2.0mm	Ø 2.7mm	Ø 2.0mm	Ø 2.5mm	Ø 2.7mm	Ø 2.7mm	Ø 3.2mm	Ø 3.5mm	Ø 3.2mm	Ø 3.5mm

SPIRAL DRILL	Ø 1.1mm	Ø 1.8mm	Ø 1.5mm	Ø 1.8mm	Ø 2.0mm	Ø 2.4mm	Ø 2.0mm	Ø 2.7mm	Ø 2.0mm	Ø 2.5mm	Ø 2.7mm	Ø 2.7mm	Ø 3.2mm	Ø 3.5mm	Ø 3.2mm	Ø 3.5mm
<b>Hand Locking Plates System</b>																
<b>Ulna Osteotomy Locking Plate II</b>																
<b>Distal Ulna Locking Plate</b>																
<b>PROlock Radius Locking Plate II</b>																
<b>PROlock Radius Locking Plate PHOENIX</b>																
<b>Tubular Locking Plates System</b>																
<b>Hallux Osteotomy Locking Plate</b>																
<b>Foot Locking Plates System</b>																
<b>Clavicle Locking Plate System NARROW</b>																
<b>Clavicle Locking Plate System</b>																
<b>Proximal Flanges Humeral Locking Plate PHOENIX</b>																
<b>Proximal Standard Humeral Locking Plate PHOENIX</b>																
<b>Distal Humeral Locking Plates</b>																
<b>Olecranon Locking Plate System</b>																
<b>Straight Compression Locking Plate</b>																
<b>Straight Locking Plate System*</b> <small>* except Straight Locking Plate 3.5/4.5mm</small>																
<b>Universal Forearm Locking Plates System</b>																
<b>Pelvic Small Fragment Reconstruction System PHOENIX</b>																
<b>PROximal MEDial Tibia Locking Plate</b>																
<b>Pilon Locking Plates System</b>																
<b>Distal Anterolateral Tibia Locking Plate</b>																
<b>Fibula Locking Plate</b>																
<b>CALcaneus Locking Plate</b>																
<b>Straight Locking Plate 3.5/4.5mm</b>																
<b>Pelvic Large Fragment Reconstruction System PHOENIX</b>																
<b>Distal Femur Locking Plate</b>																
<b>Proximal Lateral Tibia Locking Plate</b>																

SCREWDRIVER

Torque T5 

Torque T9 

Hexagon WS 2.5 

Hexagon WS 3.5 