

SMALL Fragment

**Distal
Humeral
Locking Plates**



**Straight
Compression
Locking Plate**



**CALcaneus mi
Locking Plate**



**Clavicle
Locking Plate
System**



**Olecranon
Locking Plate
System**



**Universal
Forearm Locking Plates
System**



**CALcaneus
Locking Plate**



**Anterior
Clavicle
Locking Plates
System**



**Straight
Locking Plate
System
1.5/2.0mm**



**Fibula
Locking Plate**



**Pilon
Locking Plates
System**



**Proximal
Humeral
Locking Plate
PHOENIX**



**Pelvic Small Fragment
Reconstruction
System
PHOENIX**




**Distal Anterolateral
Tibia
Locking Plate**



**PROximal
MEDial
Tibia Locking Plate**



Ø 3.5mm Locking
Cortical Screw
(Ø 2.7mm Spiral Drill)



Ø 3.5mm Locking
Cancellous Screw
(Ø 2.0mm Spiral Drill)




In Situ Benders



Screwdriver SW 2.5
or Hexagon Shank



Ø 3.5mm
Cortical Screw
(Ø 2.7mm Spiral Drill)



Ø 4.2mm Locking
Cancellous Screw
(Ø 2.5mm Spiral Drill)




Plate Holder



LARGE Fragment

**Straight
Locking Plate
3.5/4.5mm**



**Distal
Femur
Locking Plate**



**Pelvic Large Fragment
Reconstruction
System
PHOENIX**



**Proximal Lateral
Tibia
Locking Plate**



Ø 4.5mm
Locking
Cortical Screw
(Ø 3.2mm
Spiral Drill)



Ø 4.5mm
Cortical
Screw
(Ø 3.2mm
Spiral Drill)



Ø 5.9mm
Locking
Cancellous Screw
(Ø 3.5mm
Spiral Drill)



Ø 5.9mm
Cancellous
Screw
(Ø 3.5mm
Spiral Drill)



Screwdriver SW 3.5
or Shank SW 3.5



MINI Fragment

**PROlock
Radius
Locking Plate
PHOENIX**



**PROlock
Radius
Locking Plate II**



**Ulna
Osteotomy
Locking Plate II**



**Foot
Locking Plates
System**



Ø 2.7mm
Cortical Screw
(Ø 2.0mm
Spiral Drill)




Ø 3.0mm Locking
Cortical Stabilization
Screw
(Ø 2.4mm
Spiral Drill)



Ø 3.0mm Locking
Cancellous
Stabilization Screw
(Ø 2.0mm
Spiral Drill)



Ø 2.4mm Locking
Stabilization Screw
(Ø 1.8mm
Spiral Drill)



Screwdriver, Torque T9x70
or Torque-Shank




Plate Holder



MICRO Fragment

**Hand
Locking Plates
System**




Ø 1.5mm Locking
Stabilization Screw
(Ø 1.1mm Spiral Drill)



Ø 2.3mm Locking
Stabilization Screw
(Ø 1.8mm Spiral Drill)



Ø 1.8mm Cortical Screw
(Ø 1.5mm Spiral Drill)



Torque-Shank T5x80 with Handle,
AO-Connector



**Intramedullary
Nailing
System
PROXIMAL FEMUR**



Short | Intermediate | Long
180 - 480mm Length
3 CCD Angle = 120° | 125° | 130°

Lag Screw IF0-I0-xxA-S

DyCon® Set Screw IF0-7-xxA-S

Locking Screw IUI-5-xxA-S

**Flexible
Children
Nail**



Ø 3.4mm
Cortical
Screw 32346-xx

Handle with Small AO-Connector

Insertion
Extraction

Striking Rod

SCREWS

**Cannulated
Screws**

Ø 8.5mm Cannulated
Cancellous Screw
(Ø 6.2mm
Spiral Drill)



Screw Driver, SW 2.5 - 5.0
Cannulated, Shank

Guide Wire, Ø 1.2 - 3.2mm

**Cannulated
Recon
Cancellous
Screw**



Torque-Shank T25,
Handle,
ZIH-Connector



**Headless
Compression
Screws**

Ø 2.0mm/2.5mm
(L=10mm - 30mm)
Cannulated
Compression Screw



Ø 4.3mm
(L=16mm - 60mm)
Cannulated
Compression Screw




Torque-Shank or Hexagon Shank,
Handle, AO-Connector



**Twist-
Off
Screw**

Ø 2.0mm/2.7mm
Twist-Off Screw
L=11-22mm



Invers-Phillips-Screwdriver
Twist-Off Screw
AO Silicone Handle

