

ITS.

Implants for Trauma Surgery

Cannulated for minimally
invasive treatment and
guided insertion


Thread lead difference between
head and shaft thread to optimize
the fracture or osteotomy
compression

Special head design for
easier countersink of
the screw

Self-drilling &
self-tapping

Back-tapping flank for
easier explantation

Headless Compression Screws





Headless Compression Screws

Various dimensions for a wide
range of indications:

2.0, 2.5, 3.0, 3.5, 4.3, 7.5mm

- Small, space-saving trays / sets
- Intuitive / well arranged instruments



Disclaimer:

The intended users are limited to medical personnel with appropriate product training by the medical product consultants or knowledge of the surgical procedure to be applied. The medical staff must ensure that the use of I.T.S. GmbH medical devices is appropriate, taking into account the medical condition and medical history of the patient. Prior to product use, medical personnel must refer to complete information on product label and in IFU, including, but not limited to, indications, contraindications, warnings and preventative measures, and cleaning and sterilization instructions. Product availability is dependent on country registrations and clearances. For more information, please visit www.its-implant.com or contact us at office@its-implant.com. All information herein is the intellectual property of I.T.S. GmbH.

Order: No. HCS-IP-0523-US
Edition: May/2023



HEADQUARTER I.T.S. GmbH

Autal 28, 8301 Lassnitzhöhe, Austria
Tel.: +43 (0) 316/ 211 21 0
office@its-implant.com
www.its-implant.com

I.T.S. Deutschland GmbH
Thurn-und-Taxis-Platz 6, 60313 Frankfurt
Tel.: +49 69 25 73 75 322
germany@its-implant.com
www.its-implant.com

I.T.S. USA
1778 Park Avenue N - Suite 200
Tel.: 407-971-8054
info@its-implantusa.com
www.its-implant.com

ITS.
Implants for Trauma Surgery



www.its-implant.com



ITS. GmbH



ITS. Implants



[its.implants](https://www.instagram.com/its.implants)