

| | Cancellous Screw | Cancellous Screw | Cortical Screw | Stabilization Screw | Cancellous Stabilization Screw | Cortical Stabilization Screw | Cortical Screw | Cortical Screw | Cancellous Screw | Cancellous Screw | Cortical Screw | Cortical Screw | Cortical Screw | Cancellous Screw | Cortical Screw | Cancellous Screw |
|---------------------|------------------|------------------|----------------|---------------------|--------------------------------|------------------------------|----------------|-----------------|------------------|------------------|----------------|----------------|-----------------|------------------|----------------|------------------|
| | locking Ø 1.5mm | locking Ø 2.3mm | Ø 1.8mm | locking Ø 2.4mm | locking Ø 3.0mm | locking Ø 3.0mm | Ø 2.7mm | locking Ø 3.5mm | locking Ø 3.5mm | locking Ø 4.2mm | Ø 3.5mm | Ø 3.5mm | locking Ø 4.5mm | locking Ø 5.9mm | Ø 4.5mm | Ø 5.9mm |
| | 3715l-xx | 3725l-xx | 3215l-xx | 3724l-xx | 37303-xx | 37304-xx | 3227l-xx | 3735l-xx-N | 37352-xx-N | 37422-xx-N | 3235l-xx | 32352-xx | 37455-xx | 37592-xx | 32455-xx | 3059l-xx |
| SPIRAL DRILL | Ø 1.1mm | Ø 1.8mm | Ø 1.5mm | Ø 1.8mm | Ø 2.0mm | Ø 2.4mm | Ø 2.0mm | Ø 2.7mm | Ø 2.0mm | Ø 2.5mm | Ø 2.7mm | Ø 2.7mm | Ø 3.2mm | Ø 3.5mm | Ø 3.2mm | Ø 3.5mm |

| | | | | | | | | | | | | | | | | | |
|---|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
| Hand Locking Plates System | | | | | | | | | | | | | | | | | |
| Ulna Osteotomy Locking Plate II | | | | | | | | | | | | | | | | | |
| Distal Ulna Locking Plate | | | | | | | | | | | | | | | | | |
| PROlock Radius Locking Plate II | | | | | | | | | | | | | | | | | |
| Tubular Locking Plates System | | | | | | | | | | | | | | | | | |
| Hallux Osteotomy Locking Plate | | | | | | | | | | | | | | | | | |
| Foot Locking Plates System | | | | | | | | | | | | | | | | | |
| Clavicle Locking Plate System NARROW | | | | | | | | | | | | | | | | | |
| Clavicle Locking Plate System | | | | | | | | | | | | | | | | | |
| Proximal Flanges Humeral Locking Plate PHOENIX | | | | | | | | | | | | | | | | | |
| Proximal Standard Humeral Locking Plate PHOENIX | | | | | | | | | | | | | | | | | |
| Distal Humeral Locking Plates | | | | | | | | | | | | | | | | | |
| Olecranon Locking Plate System | | | | | | | | | | | | | | | | | |
| Straight Compression Locking Plate | | | | | | | | | | | | | | | | | |
| Straight Locking Plate System* <small>* except Straight Locking Plate 3.5/4.5mm</small> | | | | | | | | | | | | | | | | | |
| Universal Forearm Locking Plates System | | | | | | | | | | | | | | | | | |
| Pelvic Small Fragment Reconstruction System PHOENIX | | | | | | | | | | | | | | | | | |
| PROximal MEDial Tibia Locking Plate | | | | | | | | | | | | | | | | | |
| Pilon Locking Plates System | | | | | | | | | | | | | | | | | |
| Distal Anterolateral Tibia Locking Plate | | | | | | | | | | | | | | | | | |
| Fibula Locking Plate | | | | | | | | | | | | | | | | | |
| CAlcaneus Locking Plate | | | | | | | | | | | | | | | | | |
| Straight Locking Plate 3.5/4.5mm | | | | | | | | | | | | | | | | | |
| Pelvic Large Fragment Reconstruction System PHOENIX | | | | | | | | | | | | | | | | | |
| Distal Femur Locking Plate | | | | | | | | | | | | | | | | | |
| Proximal Lateral Tibia Locking Plate | | | | | | | | | | | | | | | | | |

| | | | | |
|--------------------|------------------|------------------|-----------------------|-----------------------|
| SCREWDRIVER | Torque T5 | Torque T9 | Hexagon WS 2.5 | Hexagon WS 3.5 |
|--------------------|------------------|------------------|-----------------------|-----------------------|