




SMALL Fragment

Distal Humeral Locking Plates 	Straight Compression Locking Plate 	CALcaneus mi Locking Plate 
Clavicle Locking Plate System 	Olecranon Locking Plate System 	Universal Forearm Locking Plates System 
Anterior Clavicle Locking Plates System 	Straight Locking Plate System 1.5/2.0mm 	CALcaneus Locking Plate 
Proximal Humeral Locking Plate PHOENIX 	Pelvic Small Fragment Reconstruction System PHOENIX 	Fibula Locking Plate 
	Distal Anterolateral Tibia Locking Plate 	Pilon Locking Plates System 
		PROximal MEDial Tibia Locking Plate 

Ø 3.5mm Locking Cortical Screw (Ø 2.7mm Spiral Drill) 	3735I-xx-N 	Ø 3.5mm Locking Cancellous Screw (Ø 2.0mm Spiral Drill) 	37352-xx-N 	In Situ Benders 	Screwdriver SW 2.5 or Hexagon Shank 
Ø 3.5mm Cortical Screw (Ø 2.7mm Spiral Drill) 	3235I-xx 32352-xx 	Ø 4.2mm Locking Cancellous Screw (Ø 2.5mm Spiral Drill) 	37422-xx-N 		Plate Holder 

LARGE Fragment

Straight Locking Plate 3.5/4.5mm 	Distal Femur Locking Plate 
Pelvic Large Fragment Reconstruction System PHOENIX 	Proximal Lateral Tibia Locking Plate 

Ø 4.5mm Locking Cortical Screw (Ø 3.2mm Spiral Drill) 	37455-xx 	Ø 4.5mm Cortical Screw (Ø 3.2mm Spiral Drill) 	32455-xx 
Ø 5.9mm Locking Cancellous Screw (Ø 3.5mm Spiral Drill) 	37592-xx 	Ø 5.9mm Cancellous Screw (Ø 3.5mm Spiral Drill) 	3059I-xx 
Screwdriver SW 3.5 or Shank SW 3.5 			

MINI Fragment

PROlock Radius Locking Plate PHOENIX 	Distal Ulna Locking Plate 
PROlock Radius Locking Plate II 	Hallux Osteotomy Locking Plate 
Ulna Osteotomy Locking Plate II 	Tubular Locking Plates System 

Foot Locking Plates System

A collection of various locking plates and screws. At the top, there are several long locking plates of different shapes, some with multiple locking holes. Below them are several individual screws of different sizes and types, including cortical screws and locking screws.

Ø 2.7mm
Cortical Screw
(Ø 2.0mm
Spiral Drill)

3227I-xx

A cortical screw with a double-flute design, used for securing plates to bone.

Ø 3.0mm Locking
Cancellous
Stabilization Screw
(Ø 2.0mm
Spiral Drill)

37303-xx

A locking cancellous stabilization screw with a wide, flat head and a threaded shaft.

Ø 3.0mm Locking
Cortical Stabilization
Screw
(Ø 2.4mm
Spiral Drill)

37304-xx

A locking cortical stabilization screw with a wide, flat head and a threaded shaft.

Ø 2.4mm Locking
Stabilization Screw
(Ø 1.8mm
Spiral Drill)

3724I-xx

A locking stabilization screw with a wide, flat head and a threaded shaft.


Screwdriver, Torque T9x70
or Torque-Shank

A black-handled screwdriver with a silver-colored shaft and a silver-colored tip.

Plate Holder


A long, thin, silver-colored metal plate holder with a small handle at one end.

MICRO Fragment


Hand Locking Plates System 
--

Ø 1.5mm Locking Stabilization Screw (Ø 1.1mm Spiral Drill) 	3715I-xx 
Ø 2.3mm Locking Stabilization Screw (Ø 1.8mm Spiral Drill) 	3725I-xx 
Ø 1.8mm Cortical Screw (Ø 1.5mm Spiral Drill) 	3215I-xx 
Torque-Shank T5x80 with Handle, AO-Connector 	

NAILING System

Intramedullary Nailing System PROXIMAL FEMUR		
Short Intermediate Long 180 - 480mm Length 3 CCD Angle = 120° 125° 130°		
Lag Screw		IF0-I0-xxA-S 
DyCon® Set Screw		IF0-7-xxA-S 
Locking Screw		IUI-5-xxA-S 

SCREWS

CAnnulated Screws 	
Ø 8.5mm Cannulated Cancellous Screw (Ø 6.2mm Spiral Drill) 	Ø 2.0mm - 7.3mm Cannulated Screw (Ø 2.6 - 5.0mm Spiral Drill) 
Screw Driver, SW 2.5 - 5.0 Cannulated, Shank 	Guide Wire, Ø 1.2 - 3.2mm 
Cannulated Recon Cancellous Screw 	3I75I-xx 
Ø 7.5mm Cannulated Recon Cancellous Screw, L=40-140mm (Spiralbohrer Ø 5.5mm/Ø 2.6mm L=180mm, Kanüliert) 	
Torque-Shank T25, Handle, ZIH-Connector 	

Headless Compression Screws 	Ø 2.0mm/2.5mm (L=10mm - 30mm) Cannulated Compression Screw 	3I257-xx 	3I207-xx 
Ø 4.3mm (L=16mm - 60mm) Cannulated Compression Screw 	3I437-xx 	Ø 3.0mm (L=10mm - 40mm) / 3.5mm (L=14mm - 30mm) Cannulated Compression Screw 	3I358-xx 
Ø 7.5mm (L=35mm - 120mm) Cannulated Compression Screw 	3I757-xx 	Torque-Shank or Hexagon Shank, Handle, AO-Connector 	
Twist-Off Screw 	Ø 2.0mm/2.7mm Twist-Off Screw L=11-22mm 	3327I-xx 	3320I-xx 
Invers-Phillips-Screwdriver Twist-Off Screw AO Silicone Handle 			