# • ITS. Values

**ITS.** stands for Implants for Trauma Surgery and represents the pioneers of a new generation. We empower surgeons and operational staff with pioneering products made in Austria, offering implants that boast full compatibility and intuitive ease of use for superior results.

### Uncompromising commitment.

All **ITS**. implants are crafted from 100% titanium, ensuring exceptional strength and biocompatibility. **ITS**. offers comprehensive sets that cover over 90% of all trauma indications—not just the most profitable areas, providing the necessary tools for achieving superior outcomes across a wide range of procedures.

### Developed with healthcare professionals.

Immediate decision processes lead to initial prototypes. After the highest clinical and biomechanical tests are performed, the final implant is released. We apply the highest global quality standards to all processes and products and renowned surgeons from around the world are involved in the development of the products.

### The Art of Trauma Surgery.

**ITS.** strives for innovation, consistently taking a different approach to trauma surgery. This distinct perspective is reflected through the artwork of Oskar Stocker, elevating our implants to the level of art.

## • Titanium • Polyaxial • Locking

All **ITS**. plates are made from Titanium Grade 2, whereas the screws are made of a harder titaniumalloy. This allows the plates to have only non-threaded holes, which all (with the exception of oblong holes) accept both locking and non-locking screws.

- Multi-directional locking (+/- I5° Locking)
- Anatomical design
- Plate material: Titanium
- Material of screws: TiAl6V4 ELI
- Easier removal of implant after fracture has healed
- ANODIZATION Type II anodized implant surfaces Improved fatigue strength of implant Reduced risk of cold welding Reduced risk of inflammation and allergy

# • Product Training

We also support our customers in the use of products and exclusively make available specialist staff who are competent with the processes and procedures of the operating theatre.

At **ITS**. we attach great importance to knowledge transfer. importance to knowledge transfer. In addition to practical training on the preparation for medical staff in our own training center at the Autal site, we also organize conferences and symposia. We invite participants to enable an intensive exchange of opinions and exchange of opinions and experiences.



**ITS**. is in business to pioneer an unparalleled user experience for surgeons worldwide. Always going beyond conventional.

# THE ART of TRAUMA SURGERY

#### Disclaimer:

Prior to product use, medical personnel must refer to complete information on product label and in IFU, including, but not limited to, indications, contraindications, warnings and preventative measures, and cleaning and sterilization instructions. Product availability is dependent on country registrations and cleannes. For more information, please visit www.its-implant.com or contact us at office@its-implant.com. All information herein is the intellectual property of I.T.S. GmbH.

HEADQUARTER

I.T.S. GmbH Autal 28, 830I Lassnitzhöhe, Austria Tel.: +43 (0) 3i6/ 2il 2i 0 office@its-implant.com www.its-implant.com

> I.T.S. USA I778 Park Avenue N - Suite 200 Maitland, FL 3275I Tel.: 407-97I-8054 info@its-implantusa.com www.its-implant.com

Order No. Z-Folder-0525-US Edition: May/2025

© I.T.S.GmbH Graz/Austria 2025 Subject to technical alterations, errors and misprints excepted.

# ITS. IMPLANTS FOR TRAUMA SURGERY



# Made in Austria



ITS, GmbH ITS, Implants its, implant

www.its-implant.com

# Upper Extremity

### • Lower Extremity

• Screws

